

Luxeo 4Z Stereo Zoom Microscope

**LUXEO 4Z
STEREO
(ZOOM)**



4144000

Stand: Dual light stable stand with rack and pinion focusing controls

Viewing Body: Binocular tube, 45°, IPD 54-74mm

Eyepieces: Wide Field eyepieces WF 10x/20mm focusable paired with eyeguards

Zoom ratio: 4.4 : 1

Magnification: 8x-35x

Illumination: Four way incident light setting with LED's for true oblique illumination, Rheostat on ring light enables cotinuous intensity control, 3-LED base light for transmitted illumination.

Electrical: Universal power supply (110V-240V, 50/60 Hz).

Packed in styrofoam box, with operation manual, allen wrench, Nylon dust cover, soft cleaning cloth, power cord and spare fuse

Objectives	WD (mm)	Eyepieces					
		W.F. 10x/20		W.F.16x/16		W.F. 20x/12	
		Field (mm)	Mag	Field (mm)	Mag	Field (mm)	Mag
none	97	6.6 to 28.5	8x to 35x	5.3 to 22.8	12.8x to 56x	4 to 17.3	16x to 70x
1.5x	45	4.4 to 19.0	12x to 52.5x	3.5 to 15.2	20.5x to 84x	2.7 to 11.5	24x to 105x